



Korean Collection for Type Cultures (KCTC), Biological Resource Center (BRC)
Korea Research Institute of Bioscience and Biotechnology (KRIBB) 181 Ipsin-gil,
Jeongeup-si, Jeollabuk-do 56212, Korea Tel: 063-570-5642, FAX: 063-570-5609

STR 결과보고서 (Result of STR matching analysis)

Cell line	Locus names								
	D5S818	D13S317	D7S820	D16S539	vWA	TH01	AM	TPOX	CSF1PO
Reference (TCCSUP)	12,12	11,14	8,9	9,11	14,16	6,9.3	X,X	8,8	10,10
Your cell (TCCSUP1)	12,12	11,14	8,9	9,11	14,16	6,9.3	X,X	8,8	10,10
Note	○	○	○	○	○	○	○	○	○

STR Kit : AmpFISTR Identifier (Applied Biosystems, USA)

Analysis Methods : 3130XL DNA analyzer (Applied Biosystems, Life Technologies)
GeneMapper ID v3.2 (Applied Biosystems, Foster, CA)

Analysis Institution : HumanPass INC.

The Fingerprinting of cell lines by 'AmpFLSTR identifier PCR Amplification kit' was tested on July, 2016.

의뢰하신 cell line은 Reference 결과와 일치합니다. (Your cell matches with the reference.)



생물자원센터장
Director of KCTC



Korean Collection for Type Cultures (KCTC), Biological Resource Center (BRC)
Korea Research Institute of Bioscience and Biotechnology (KRIBB) 181 Ipsin-gil,
Jeongeup-si, Jeollabuk-do 56212, Korea Tel: 063-570-5642, FAX: 063-570-5609

STR 결과보고서 (Result of STR matching analysis)

Cell line	Locus names								
	D5S818	D13S317	D7S820	D16S539	vWA	TH01	AM	TPOX	CSF1PO
Reference (TCCSUP)	12,12	11,14	8,9	9,11	14,16	6,9.3	X,X	8,8	10,10
Your cell (TCCSUP2)	11,11	11,12	8,12	12,12	15,17	6,8	X,X	8,11	11,11
Note	X	X	X	X	X	X	O	X	X

STR Kit : AmpFISTR Identifier (Applied Biosystems, USA)

Analysis Methods : 3130XL DNA analyzer (Applied Biosystems, Life Technologies)
GeneMapper ID v3.2 (Applied Biosystems, Foster, CA)

Analysis Institution : HumanPass INC.

The Fingerprinting of cell lines by 'AmplFLSTR identifier PCR Amplification kit' was tested on July, 2016.

의뢰하신 cell line은 Reference 결과와 일치하지 않습니다. (Your cell does not match with the reference.)



생물자원센터장
Director of KCTC